

SF LIPO MED

SLIMFIT AESTHETICS

SF Lipo Med

The SF Lipo Med is a fabulous innovation, which incorporates four incredible technologies in one hand piece, which provides fantastic instant results, in the fraction of the time.

These technologies used are:

- Cavitation (Breaks down Fat)
- Radio Frequency (Tightens the Skin)
- LED (Rejuvenates the Skin)
- Vacuum (Kick-starts the Lymphatic System)

With its high-powered intensity reaching an incredible 100 KHz, this device is likely to work wonders, as it is one of the most effective machines available to date, and can be used all over the body and face. The treatments are excellent for fat reduction, skin tightening, and aid in reducing stretch marks and cellulite. Within one 20-minute session, the fat cells are broken down and the skin tightened, resulting in up to 4 inches of fat loss.

This is the perfect solution for any organisation that is looking to provide non-surgical alternatives and expand their clientele.




www.slimfitaesthetics.com

[f /slimfitaesthetics](https://www.facebook.com/slimfitaesthetics)

[i /slimfitaesthetics](https://www.instagram.com/slimfitaesthetics)

[g+ /slimfitaesthetics](https://www.google.com/+slimfitaesthetics)



How it Works

Cavitation (Breaking down Fat):

- Cavitation uses ultrasound (positive and negative alternation of sound waves) to disrupt the fat cell membranes. Usually a fat cell membrane cannot withstand this pressure, it will implode, and turn into fatty acids.
- The fatty liquid is used as an energy source and the rest of the insoluble fatty acids are transported to the liver and are processed.

Radio Frequency (Tightening the Skin):

- Radio Frequency heats the deeper dermis, to a constant 40-45 degrees, stimulates fibroblast action, which are connective tissues that produce new collagen and elastin fibres.
- This causes an immediate contraction of the fibres, therefore the skin is instantly tighter.

LED (Rejuvenating the Skin):

- LED lights speeds cell production, making skin regenerate faster, which will result in more youthful, and fresher-looking skin.

Vacuum (Kick start the Lymphatic System):

- Vacuum suction treatment is applied to the surface of the skin tissue causing lift and stimulation of the under-lying structures.
- Vacuum suction is aimed at improving blood flow and kick-starting the lymphatic system, thus flushing out toxins through the lymph nodes.
- The vacuum opens the lymph nodes to provide lymphatic drainage that helps the body to eliminate the waste products that have been produced during the treatment.

ALL ABOVE TECHNOLOGIES ARE COMBINED IN JUST ONE SPECIALLY-DESIGNED HANDPIECE!

Key Attributes:

- Instant Inch Loss
- Cellulite Reduction
- Skin Tightening
- Stretch Marks Reduction
- Skin Rejuvenation
- Wrinkle Reduction
- Facial Lifting
- Breast Lifting
- Brazilian Bum Lifting
- Increased Metabolism
- Improved Blood Circulation



SF LIPOMED
SLIMFIT AESTHETICS

www.slimfitaesthetics.com

 /slimfitaesthetics

 /slimfitaesthetics

 /slimfitaesthetics